

Other NGS applications: CNVs, Structural Variations

Erika Souche

Definitions

- Structural Variant (SV)


- Copy Number Variant (CNV)

- Deletion
 - Duplication
 - Triplication
 - ...


- Insertion

- Inversion


- Translocation



CNV signatures




CNV signatures




Detecting CNVs

- Based on Read Counts (RC)




Detecting CNVs

- Based on Read Counts (RC)




Detecting CNVs

- Based on Read Counts (RC)




Detecting CNVs

- Based on Read Counts (RC)




Detecting CNVs

- Based on Read Counts (RC)




Detecting CNVs

- Based on Read Counts (RC)




Detecting CNVs

- Based on Read Pairs (RP)




Detecting CNVs

- Based on Split Reads (SR)




Detecting CNVs

- Based on Split Reads (SR)




Detecting CNVs

- Deletion supported by
 - Read counts
 - Read pairs
 - Split reads



Detecting SVs



Summary

- SV detection possible by WGS
- Resolution dependent on
 - Depth of coverage
 - Detection algorithm
- Limitations
 - Eventual need of reference
 - SVs in segmental duplications, repetitive regions
 - Complex SVs

SVs in targeted sequencing?

- Strategies

- Read Counts
- Read Pairs
- Split Reads


} Requires reference set
} Only if breakpoints are part of/close to targeted regions

- Limitations

- Bias introduced during capture/PCR
- Resolution dependent on targeted regions
- No information on breakpoints


SVs in targeted sequencing

- Capture sequencing



SVs in targeted sequencing

- Capture sequencing



SVs in targeted sequencing

- Capture sequencing


- RPKM

- Reads per thousand bases per million reads sequenced


- 1 RPKM for each capture target/probe

$$RPKM = \frac{Read\ starts * 10^9}{Total\ reads * Target\ length}$$

SVs in capture sequencing




SVs in capture sequencing




SVs in capture sequencing


RPKM by probe



SVs in capture sequencing



SVs in capture sequencing



SVs in targeted sequencing


- Amplicon sequencing
 - RPKM
 - Deviations detected by
 - Z scores
 - Log2 ratios

$$Z = \frac{RPKM - \text{Mean}}{\text{Standard deviation}}$$

$$\log_2 = \log\left(\frac{RPKM}{\text{Mean}}\right) / \log(2)$$


	Deletion	Duplication
Z score	< -3	> 3
Log2 ratio	< -0.7	> 0.5

SVs in amplicon sequencing




SVs in amplicon sequencing

Gene	Transcript	Exon	Amplicons/exon	Average log ratio	Average Z score
LDLR	NM 000527.4	7	2	-1.0713	-13.4763
		8	2	-1.0569	-5.3929



SV filtering

- Long lists of variants
 - Filtering required to discard False Negatives
 - Read ratio
 - CNV scores
 - CNV length
 - Etc.



Summary

- SVs can be called from short read data
- Some SVs are more difficult/impossible to detect from targeted data
- Other NGS strategies
 - *de novo* assembly
 - Long read sequencing

Questions?

Erika Souche